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6. Annexe : T. G. Appleton *J. Chem. Educ.* **1977**, *54*, 443.

Trevor G. Appleton
University of Queensland
Brisbane 4067
Australia

Oxygen Uptake by a Cobalt(II) Complex

An undergraduate experiment

Natural oxygen-carrier and storage proteins contain a transition metal to which the oxygen reversibly coordinates—e.g., iron (myoglobin, hemoglobin), copper (hemocyanin), vanadium (hemovanadin). Various simpler coordination compounds will also reversibly coordinate oxygen, and have been extensively studied as “model compounds” (1, 2). Cobalt(II) complexes are among the most studied. Two types of bonding of dioxygen to cobalt(II) have been distinguished





Figure 1. Structure of the oxygen adduct.

That is, the Co:O₂ ratio may be 1:1 or 2:1, depending on the nature of the ligands, temperature, solvent, etc. (2-7).

One of the best known cobalt(II) “oxygen carriers” is Co(salen) (I), where salenH₂ represents N,N'-bis(salicylaldehyde)ethylenediimine

Co(salen) exists in two different solid forms, depending on method of preparation: a brown form *active* toward absorption of oxygen, and a dark red *inactive* form. The *inactive* form consists of dimeric units (II) in which the Co atom of one Co(salen) molecule interacts with an O-atom of the second (8). The *active* form also has dimeric units, but with one Co atom directly above the other (III) (9). Channels in the structure allow oxygen to enter easily. The *active* form absorbs oxygen at room temperature, and releases it at higher temperatures. This cycle may be repeated many times, although the activity of the compound toward oxygen uptake slowly decreases on continued cycling, owing to decomposition (10).

When Co(salen) is dissolved in donor solvents (e.g., dimethylsulfoxide DMSO, dimethylformamide DMF, pyridine) in the presence of oxygen, adducts rapidly form (6, 7). In DMSO, a dark brown solid, [(DMSO)Co(salen)]₂O₂ precipitates, which probably has a structure similar to that of the DMF analog (Fig. 1, L - DMF) (11). The dia-


Figure 2. Preparation of Co(salen).



Figure 3. Oxygen uptake by Co(salen).

magnetic compounds may be formally regarded as containing two Co(III) atoms linked by a peroxo bridge.

In this experiment, the inactive form of Co(salen) is prepared by an adaptation of the method of Bailes and Calvin (5). An aqueous solution of cobalt(II) acetate is added to a hot ethanol solution of salenH₂ under nitrogen (Fig. 2). A mixture of active and inactive forms is initially obtained, which is entirely converted to the inactive form when heated in ethanol. Once inactive Co(salen) has been formed, it may be filtered off and handled in air.

Oxygen uptake is measured using the simple apparatus depicted in Figure 3. The side-arm test tube is flushed with oxygen to increase the rate of reaction, and to leave no doubt that the gas absorbed is, in fact, oxygen. It is not necessary to remove all traces of nitrogen. When the side-arm tube is inverted, allowing the DMSO to dissolve the Co(salen), oxygen is rapidly absorbed and a dark brown precipitate forms. Reaction is complete in 10–20 min, if the tube is continuously shaken. The small volume of the apparatus prevents small temperature changes from having a large effect on volume. In the two years that this experiment has been used in a third-year Bioinorganic Chemistry course at the University of Queensland, most student results for the O₂:Co ratio have been close to the theoretical value, 0.5, and all have been within the limits 0.44–0.60.

The oxygen adduct formed is a very fine precipitate, which cannot be easily filtered, but may be isolated by centrifuging. When chloroform is added to the adduct, it slowly dissolves, releasing fine streams of oxygen bubbles, to give a red solution of Co(salen). The student can be asked to discuss the factors which affect the equilibrium (3) in DMSO and chloroform solvents

Electron donation by a ligand, *L*, is required to stabilize the structure in Figure 1. Chloroform is a poor donor unlikely to coordinate to cobalt. In DMSO, the high concentration of DMSO provided by the solvent, and the low solubility of the adduct help to force the equilibrium to the right.

The total time required for experimental work, including characterization of compounds by IR and melting/decomposition behavior,

is about 6 hr, easily broken if necessary into two 3-hr sessions.

Experimental

Preparation of *N,N'*-Bis(salicylaldehyde)ethylenediimine, salenH₂

To a solution of 4.9 g (4.2 ml) salicylaldehyde in 50 ml boiling 95% ethanol add 1.2 g (1.35 ml) ethylenediimine. Stir the reaction mixture for 3½ min, and leave the solution to cool in an ice bath. Filter the bright yellow flaky crystals under suction, wash with a small volume of ethanol, and air-dry. Melting point 120°C. Yield 90–100%.

Preparation of the Inactive Form of (*N,N'*-bis(salicylaldehyde)ethylenediimine)cobalt(II), Co(salen)

Weigh 2.34 g salenH₂ into a 250-ml 3-necked flask. Add 120 ml 95% ethanol. Set up the apparatus as shown in Figure 2. Stir using a magnetic follower, and flush the apparatus with nitrogen. Adjust the nitrogen flow to a steady rate (~one bubble/5 s). Immerse the flask in a water bath maintained at 70–80°C. Dissolve 2.17 g cobalt(II) acetate, Co(CH₃CO₂)₂·4H₂O in 15 ml hot water. When the salenH₂ has all dissolved, add the cobalt acetate solution through the funnel. A brown gelatinous precipitate forms immediately. Continue heating and stirring for a further hour, during which time the precipitate turns dark red. Cool the flask by immersing in cold water. Discontinue the nitrogen flow and filter off the solid on a sintered glass funnel (in the air). Wash three times with 5 ml water, then with 5 ml 95% ethanol. Dry on the funnel, then in a vacuum desiccator. Decomp. ~300°C. Yield 70–90%.

Oxygen Uptake by Co(salen)

Weigh out between 0.05 and 0.1 g ground Co(salen) into a side-arm test tube (2 × 18 cm). Place approximately 5 ml DMSO in a beaker, and bubble oxygen through for a few seconds, to saturate it with oxygen (caution—DMSO is absorbed through the skin). Transfer the DMSO into a small test tube (1 × 7.5 cm) until it is filled to about 2 cm from the rim. Using forceps, carefully lower the small test tube into the side-arm tube without spillage.

Connect a length of glass tubing (graduated 0–10 ml—a pipet would suffice although the apparatus shown is easily made), to a second glass tube ("movable arm") using tubing. Arrange the tubes in a U-shape (Fig. 3) and fill with water. Connect the graduated tube to the side-arm tube with polyethylene tubing. Adjust the movable arm so that the water level is near the bottom of the graduated tube.

Flush the side-arm tube with a gentle stream of oxygen. Insert a tight-fitting rubber stopper in the mouth of the tube. Adjust the movable arm to make the water levels the same in both tubes (i.e. pressure within the apparatus is atmospheric). Note the water level in the graduated tube. Carefully invert the side-arm tube (holding near the stopper to minimize heating by the hand) so that the DMSO is introduced onto the Co(salen), but not spilled into the polyethylene tube. Gently shake the tube. As oxygen is absorbed the water level in the graduated tube begins to rise. Note the changes occurring in the tube. Continue shaking until no further change in water level occurs (10–20 min). Adjust the movable arm so that the water levels in the two tubes are again equal, and read off the new level in the graduated tube.

From the decrease in volume at room temperature and atmospheric pressure, calculate the number of moles of dioxygen absorbed per mole of Co(salen).

Reaction of the Oxygen Adduct in Chloroform

Remove the stopper from the side-arm tube, and remove as much as possible of the dark brown suspension into a centrifuge tube. Centrifuge until the precipitate has settled to the bottom of the tube (15–30 min). Carefully remove the supernatant DMSO. To the residue in the tube (drying is not necessary) add 5–10 ml chloroform without stirring. Observe the result, and attempt to explain your observations.

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6. Annexes

6.1. J.D. Woolins « Inorganic Experiments »

Dissolve 0.5 g of the iron EDTA complex in 20 cm³ nitrogen degassed water in a flask equipped with a magnetic stirrer bar and a nitrogen inlet. Add 0.5 g sodium azide and stir the solution at about 70 °C for 1 hour whilst maintaining a steady flow of nitrogen through the solution. Reduce the volume to about 10 cm³. This can be done under reduced pressure or by evaporation under a flow of nitrogen. Cool to room temperature and add 50 cm³ of ethanol. If the product precipitates as an oil, decant the supernatant and stir the residual oil with more ethanol to induce crystallisation. Record the weight obtained and calculate the percentage yield. Record the infra-red spectrum as either a KBr disc or as a Nujol mull and identify the bands due to coordinated N₂ and the CO groups of the EDTA.

Exercises

1. What is the structure of the [Fe(HY)H₂O]·H₂O complex?
2. Discuss the bonding of dinitrogen to metals.

References

- R. A. Henderson, *Adv. Inorg. Chem. Radio Chem.* **1983**, *27*, 198.
M. Hidai, *J. Amer. Chem. Soc.* **1983**, *105*, 1680.
R. J. H. Clark, *J. Chem. Soc. Dalton Trans.* **1983**, 2005.
M. Garcia Basallote, J. M. Lopez Alcalá, M. C. Puerta Vicuña, F. Gonzalez Vilchez, *Inorg. Synth.* **2007**.

3.11 Nickel Dihalide Phosphine Complexes

Ivan P. Parkin
Adapted from the Open University CHEM777 Summer School as prepared by the S343 course team

Special Safety Precautions

1. Nickel salts are toxic. They are also irritants. Wear gloves to handle all nickel salts.
2. The carbon disulphide generated in the experimental is highly toxic and flammable. This preparation must be carried out in a fume cupboard.
3. Phosphines are irritants and toxic.
4. Carry out all reactions in a fume cupboard.

The physical and chemical properties of a coordination complex are dependent on its geometry. This distinction becomes important when different geometric isomers are possible. For example, *cis* and *trans* isomers and square planar or tetrahedral isomers can have different IR spectra, magnetic moment values and UV spectra.

In this experiment, you will synthesise three nickel halide bis-phosphine complexes and determine their geometry using spectroscopic methods. From this data, you will be in a position to predict the geometry of the complexes and suggest reasons for a specific geometric isomer.

Experimental

a) NiCl₂(PPh₃)₂

Triphenylphosphine (2.8 g), propan-2-ol (30 cm³, dried) and two anti-bumping granules are placed in a dry 100 ml B24 neck round-bottomed flask. The flask is fitted with a B24 reflux condenser and the mixture brought to reflux.

Nickel dichloride (1.2 g NiCl₂·6H₂O) is dissolved in ethanol (15 ml) in a 50 ml conical flask and the solution warmed to approximately 40 °C. The warm solution of nickel dichloride is added with care down the condenser to the refluxing triphenylphosphine solution. Normally, the precipitate of NiCl₂(PPh₃)₂ will form immediately upon addition. The reaction mixture is allowed to cool to room temperature and the product filtered off using a sintered glass filter and Buchner funnel. Wash the precipitate on the filter with cold ethanol (15 ml) and cold dry ether (15 ml). Dry by suction.

Record the yield, melting point, IR (4000–200 cm⁻¹), UV (in toluene or dry CHCl₃) and magnetic moment.

b) $[\text{Ni}(\text{NCS})_4(\text{PPh}_3)_2]$

Change a B24 100 ml "Quickfit" round-bottomed flask with nickel nitrate, 1.5 g $\text{Ni}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$, ethanol (25 ml) and two antihumping granules or glass beads. Stir or swirl to dissolve, heating gently if necessary. Then add the *finely ground* sodium thiocyanate (0.8 g), top the flask with a B24 reflux condenser and reflux for 20 minutes (during which time you should prepare a phosphine solution as described in the first paragraph of the $[\text{NiCl}_2(\text{PPh}_3)_2]$ preparation). Cool the thiocyanate solution in an ice bath, while scratching the inside of the flask with a glass rod to precipitate the sodium nitrate and any unreacted sodium thiocyanate.

Filter the nickel thiocyanate solution through a sintered glass funnel into a Buchner flask. Always remove the tubing from your Buchner flask *before* you turn off the water tap, so that water does not suck back. The thiocyanate solution should be transferred to a conical flask and heated on a hotplate with one or two (not more) antihumping granules. The hot nickel thiocyanate solution is then carefully added down the reflux condenser to the refluxing phosphine solution. The reaction mixture is allowed to cool to room temperature and the product filtered off using a sintered glass filter and Buchner funnel. Wash the precipitate on the filter with cold ethanol (15 ml) and cold dry ether (15 ml). Dry by suction.

Record the yield, melting point, IR (4000–200 cm^{-1}), UV (in toluene or dry CHCl_3) and magnetic moment.

c) $[\text{NiCl}_2(\text{PCy}_3)_2]$

A dinitrogen atmosphere is essential to avoid oxidation of the phosphine. Dinitrogen is bubbled through the solution of the adduct $\text{PCy}_3 \cdot \text{CS}_2$ to remove the carbon disulphide, which is highly flammable and toxic. *This preparation must therefore be carried out in a fume cupboard.*


Figure 3.11-1. Reflux apparatus for the preparation of $[\text{NiX}_2(\text{PR}_3)_2]$.

Assemble the reflux apparatus as shown in Figure 3.11-1. There should be a small amount of liquid paraffin in the dinitrogen trap to prevent air entering the apparatus, and to allow dinitrogen to escape from the apparatus.

Flush the system with dinitrogen by passing a fairly rapid stream of dinitrogen through the apparatus for a few minutes. The pressure should not exceed 25 kPa, which gives a flow rate of about 4 bubbles per second.

Weigh 1.9 g of the tricyclohexylphosphine-carbon disulphide adduct into the flask and add dry propane-2-ol (20 ml) together with antihumping granules. Attach the condenser and pass dinitrogen through the solution for five minutes. Make sure that the dinitrogen inlet is well below the liquid surface.

Maintain the dinitrogen supply and heat the flask to reflux. If the flow is sufficient, the CS_2 will be driven off within fifteen minutes. The solution will then be clear and colourless (or a very pale straw colour), but you should continue until you are certain that there is no further colour change, and until you are ready to add the warm nickel salt. (You may have to top up with degassed solvent, since this evaporates during the reflux.)

Nickel dichloride (0.6 g $\text{NiCl}_2 \cdot 6\text{H}_2\text{O}$) is dissolved in dried ethanol (15 ml) in a 50 ml conical flask and the solution warmed to 40 °C. The warm nickel dichloride solution is added to the refluxing phosphine carefully down the condenser and the product obtained as previously described for $[\text{NiCl}_2(\text{PPh}_3)_2]$.

Record the yield, melting point, IR (4000–200 cm^{-1}), UV (in toluene or dry CHCl_3) and magnetic moment.

Exercises

- Four coordinate nickel complexes can be either square planar or tetrahedral. Draw three diagrams for $[\text{NiCl}_2(\text{PPh}_3)_2]$ showing the possible geometric isomers.
- Metal chloride stretches occur in the infra-red between 400–250 cm^{-1} . How many Ni–Cl stretches are seen for $[\text{NiCl}_2(\text{PPh}_3)_2]$ and $[\text{NiCl}_2(\text{PCy}_3)_2]$? Does this suggest a particular geometric isomer for each complex?
- Calculate the extinction coefficient ϵ for the bands in the UV spectra of the three nickel complexes. Does this value suggest tetrahedral or square planar coordination?
- Draw a simplified splitting diagram for the *d*-orbitals in a tetrahedral and square planar crystal field. Fill in the appropriate number of electrons for Ni^{2+} in both cases. Compare the number of unpaired electrons with the magnetic moments you measured. Does this predict a certain geometry?
- Compare the ligand field properties of the thiocyanate and chloride ions and use this to explain the differences in the experimental results for $[\text{NiCl}_2(\text{PPh}_3)_2]$ and $[\text{Ni}(\text{NCS})_4(\text{PPh}_3)_2]$. What are the steric and electronic differences between PPh_3 and PCy_3 , which account for your findings?
- Can you suggest any further methods we could use to differentiate between square planar and tetrahedral geometries?

3.12 Influence of Ligand Field Tetragonality on the Ground-State Spin of Nickel(II) Complexes

David M. L. Goodgame

Special Safety Precautions

N,N-Diethylethylenediamine should be stored and used in a fume cupboard and any skin contact with it avoided; wear rubber gloves when using it. Also wear rubber gloves when working with solutions of the nickel salts and their complexes, to avoid the possibility of heavy metal skin allergy.

Regularly octahedral, first series transition-metal complexes of electronic configuration d^n ($n = 4-7$) may be either high- or low-spin depending, primarily, on the strength, Δ or $10Dq$, of the ligand field. For distorted ligand fields, a wider range of possibilities exists for changes in the ground-state spin. Changes in spin can have important effects, e.g. high-spin octahedral iron(II) complexes are much more labile than their low-spin counterparts. Moreover, the haem complexes involved in haemoglobin, catalase, cytochrome C, etc., contain the iron atom in either a high- or low-spin state, depending on the axial ligand.

The effects of changes in the axial ligands, X, on the electronic properties of a tetragonally distorted complex $trans-ML_4X_2$ are well illustrated by the magnetic properties and $d-d$ spectra of a series of nickel(II) complexes of formula $Ni(Etgen)_2X_2$, where $Etgen = N,N$ -diethylethylenediamine ($Et_2NCH_2CH_2NH_2$) and $X = Cl^-$, Br^- , I^- , NCS^- , etc.

In the following experiment, a number of compounds of this type can be made and their spin-states determined by magnetic susceptibility measurements, supplemented, if possible, by their solid state (reflectance) $d-d$ spectra.

Experimental

a) $Ni(Etgen)_2(NCS)_2$

Prepare an ethanol solution of nickel thiocyanate by dissolving nickel nitrate hexahydrate (2.0 g) and finely powdered potassium thiocyanate (1.3 g) separately in hot absolute ethanol and mixing the two solutions. Keep the volume of solution to a minimum and cool the mixture thoroughly before filtering off the precipitated potassium nitrate. (For your guidance, the solubility of nickel nitrate in hot ethanol is ca. 50 g/100 ml, and that of powdered potassium thiocyanate is ca. 10 g/100 ml).

To the resulting solution of nickel thiocyanate, add 2 ml of $Etgen$ by means of a syringe, and stir or shake the solution during the addition. Filter off the precipitated complex, wash it with a little cold ethanol, and dry it (air dry or vacuum desiccator). Record the yield of crude product and then recrystallise it from methanol. Record the yield of pure, crystalline product and take its m.p.

b) $Ni(Etgen)_2I_2$

Prepare an ethanol solution of nickel iodide containing ca. 2 g NiI_2 from ethanol solutions of the appropriate amounts of nickel nitrate hexahydrate and sodium iodide (solubility of NiI_2 in hot ethanol ca. 16 g/100 ml). (Again, cool the solution well before filtering off the precipitated sodium nitrate). Add $Etgen$ (2 ml) using a syringe and filter off, wash (ethanol) and dry the precipitated complex. Do not attempt to recrystallise it. (Alternatively, or additionally, the analogous nickel nitrate complex can be made by simple addition of $Etgen$ to a concentrated ethanol solution of nickel nitrate).

c) Physical measurements

Measure the magnetic susceptibilities of the complexes and, in the case of the thiocyanate complex, also calculate the magnetic moment. If the thiocyanate complex has formed large crystals, it may be necessary to powder some of the sample to achieve good packing in the Gouy (or Evans) tube used for the magnetic measurements.

Obtain the IR spectrum of the thiocyanate complex and identify the bands due to the NCS group, which coordinates to the nickel via the nitrogen atom, ν , ρ . It is an isothiocyanato complex. If you have access to a spectrometer capable of measuring the reflectance spectra of solid samples in the UV-vis region, record the spectra of the complexes and relate the spectra to the colours of the compounds.

d) Other complexes

The analogous complexes with nickel chloride and nickel bromide exist as dihydrates as well as anhydrous compounds. Investigate the formation of these compounds by dissolving 1 g of nickel chloride or bromide respectively, in ethanol and adding $Etgen$ (1 ml) to each of these solutions. Observe the colours of the initial precipitates and any colour changes that take place when the products are filtered off, washed with ethanol and air dried. Put a small amount of each compound in a melting-point tube, record the melting point or decomposition point and note any colour changes that occur prior to that point.

If you have isolated both anhydrous and hydrated forms of each complex, record their infrared spectra (using d_2O Nujol) and identify additional bands from the water molecules. If time permits, the weight loss when one of the dihydrates is converted to its anhydrous analogue can be determined, for example, by heating in an oven or a drying pistol at ca. 100°C, or on prolonged storage in an evacuated desiccator over a good drying agent.

Exercises


Draw the crystal field splitting diagram expected for a tetragonally elongated octahedron and use it to explain why some of the compounds you have made are diamagnetic and the others are paramagnetic with magnetic moments indicative of two unpaired electrons per nickel atom.

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



6.4. R. J. Angelici « Synthesis and Technique in Inorganic Chemistry »

in characterizing a complex. Thus, the oxidation state of a metal in a complex may be readily assigned on the basis of its number of unpaired electrons. For example, copper in a given coordination compound will have 0 or 1 unpaired electron, depending on whether the oxidation state is +1 or +2. Also, the number of unpaired electrons may be helpful in assigning geometries to complexes. The geometries of 4-coordinate Ni complexes, $Ni^{II}L_4$, have been extensively studied by this method. Complexes of this formula may have either square planar or tetrahedral geometry. The d orbitals on Ni split, according to crystal field theory, in the following manner depending on the geometry of the complex:

The placement of the eight d electrons of Ni^{II} into the lowest energy orbitals leaves two unpaired electrons for tetrahedral complexes and no unpaired electrons in square planar complexes. Four-coordinate Ni^{II} complexes have been found by magnetic susceptibility studies to assume either the tetrahedral or the square planar configuration, depending on the ligands bound to the metal. There are even examples of Ni^{II} complexes that in solution, exist partially in the tetrahedral and partially in the square planar form.

Finally, the number of unpaired electrons is frequently of value in indicating the presence of unusual modes of bonding. The metal carbonyl complex, $Fe_3(CO)_9$, has long been known from x-ray structural studies to have the geometry

The coordination around each Fe atom is approximately octahedral. On the basis of this structure, one might hope to predict the magnetic

Experiment 5

Magnetic Susceptibility

Note: $[Cr(NH_3)_4(NO_2)_2]$, prepared in Experiment 4, may be used in this experiment. The procedures given in this experiment may also be used in the determination of the magnetic susceptibility of any other paramagnetic or diamagnetic substance.

The widespread use of the valence-bond theory of chemical bonding in transition metal complexes stemmed largely from its ability to account for and predict magnetic properties of a wide variety of coordination compounds. Although it has now been largely displaced by ligand field and molecular orbital theories, it provided a strong impetus for more extensive and detailed investigations of these properties. This theory accounts for the magnetism of complexes in terms of the number of unpaired electrons and their associated magnetic or spin properties. Those complexes that contain unpaired electrons are attracted into a magnetic field and are said to be paramagnetic, while those with no unpaired electrons are repelled by such a field and are called diamagnetic. Diamagnetism results from the fact that electrons in any material will move so as to generate a magnetic field that will oppose an applied field. Thus, any metal complex will possess some diamagnetism, although the overall character of the compound may be paramagnetic because of the very large paramagnetic susceptibility of an unpaired electron.

Thus, the magnitude of the interaction of a coordination compound with a magnetic field frequently allows the determination of the number of unpaired electrons in the complex. While this is usually true for complexes of the first-row transition elements, it is not always true for complexes of the second- and third-row metals. In these cases, the observed magnetic properties may reflect not only the spin properties of the electron but also that associated with the orbital motion of the electron. In this experiment, we will be studying the magnetic susceptibility of a coordination compound of a first-row transition metal, and from the result we will calculate the number of unpaired electrons in the complex.

A knowledge of its number of unpaired electrons is very helpful

properties of the compound. If the ligands are formally said to be neutral, then the Fe is in the zero oxidation state and will have eight electrons. Each Fe will bond with two electrons from each of the terminal CO groups and one electron from each of the bridging CO groups. Thus, each Fe atom is surrounded by a total of 17 bonding and valence electrons. The odd number of electrons predicts that the entire molecule would have two unpaired electrons, one on each Fe. Since the compound is, in fact, diamagnetic, a reasonable bonding proposal requires that these two electrons become paired with each other. It has therefore been suggested that there is an electron-pair, covalent bond formed between the Fe atoms (dotted line in drawing) that accounts for the diamagnetism of the complex. On re-examining the structure, the relatively short Fe-Fe distance (2.46 Å) supports this suggestion. Magnetic studies have contributed greatly to our present understanding of chemical bonding in transition metal complexes.

MAGNETIC SUSCEPTIBILITY

The purpose of this experiment is to acquaint the student with a common experimental method of measuring magnetic susceptibilities. Unfortunately, there is not space in this book to provide a background in magnetism, which should be familiar to the student. It is suggested that a physics textbook be reviewed for this information. The approach here will be to make use of equations relating to magnetic fields without attempting to justify them. A rigorous treatment may be obtained in the references given at the end of this experiment.

The magnetic field inside a substance that is placed in an external magnetic field will depend not only on the magnitude of the applied field but also upon the ability of that substance to produce its own field, which will add to (paramagnetism) or subtract from (diamagnetism) the applied field. This magnetic field, B , inside the substance may then be represented as the following sum:

$$B = H_0 + 4\pi I \quad (1)$$

where H_0 is the applied field and I is the intensity of magnetization of the substance or its induced magnetic moment. The ratio I/H_0 is a measure of the susceptibility of the substance to interact with an applied field, and is called the volume susceptibility, κ . This volume susceptibility is usually converted to susceptibility per gram of substance by dividing by its density, d . It is called the specific susceptibility, χ :

$$\chi = \kappa/d. \quad (2)$$

The magnetic susceptibility per mole of substance, χ_M , is of most value for chemical applications. It is obtained by multiplying χ by the molecular weight, M , of the compound:

$$\chi_M = (\chi)(M) \quad (3)$$

The molar susceptibility, χ_M , is positive if the substance is paramagnetic and negative if diamagnetic.

Since the overall χ_M of a compound is the sum of the susceptibilities of the paramagnetic electrons in the complex and of the diamagnetic paired electrons on the metal, ligands, and ions of the compound, the susceptibility of only the unpaired electrons, χ_M' , may be obtained from the additive relationship

$$\chi_M = \chi_M' + \chi_M(\text{metal core electrons}) + \chi_M(\text{ligands}) + \chi_M(\text{ions}) \quad (4)$$

In this equation, the last three terms represent the susceptibilities of the diamagnetic paired electrons of the inner core of the metal, the ligands, and the other ions (e.g., NO_3^- in this experiment), respectively. Thus, in order to determine the magnetic susceptibility of the unpaired electrons in the complex, it is necessary to correct the measured χ_M for the diamagnetism of the other groups. It is known that susceptibilities of diamagnetic groups change very little with their environment; hence, it is possible to calculate the diamagnetism of a molecule by adding together the diamagnetisms of structural components of the molecule. Some diamagnetic corrections for ligands and ions are given in Table 5-1. The corrections for transition metal ions have only been estimated (Table 5-1) and frequently only the sum $\chi_M' + \chi_M$ (metal core electrons) is evaluated.

TABLE 5-1. Diamagnetic Corrections for Ligands and Ions*

| Formula | $-10^6 \chi_M', \text{ cgs}$ | Formula | $-10^6 \chi_M', \text{ cgs}$ |
|----------------------------------|------------------------------|----------------------|------------------------------|
| Ag^+ | 28 | ClO_4^- | 32 |
| BF_4^- | 37 | F | 9 |
| Ba^{2+} | 23 | H_2O | 13 |
| Br | 35 | I^- | 51 |
| Ca^{2+} | 10 | K | 15 |
| CN | 13 | Mg^{2+} | 5 |
| NCS^- | 31 | NH_3 | 18 |
| CO | 10 | NO_2 | 19 |
| CO_3^{2-} | 28 | Na | 7 |
| $\text{C}_6\text{H}_6\text{O}_2$ | 30 | OH | 12 |
| $\text{C}_6\text{H}_5\text{N}$ | 49 | PH_3^+ | 32 |
| C_6H_6 | 55 | SO_4^{2-} | 40 |
| Cl | 23 | Zn^{2+} | 15 |

*The calculated inner core diamagnetism of the first-row transition metals is approximately $-13 \times 10^{-6} \text{ cgs}$.

generates a magnetic field of at least 5000 gauss. Some provision must be made to remove the magnetic field from the sample, so that the sample may be weighed in the absence and presence of the magnetic field. Since an electromagnet is commonly used, the field may be simply turned off. With a permanent magnet, the magnet must be physically moved from the sample area, in order to "turn the field off."

In determining the magnetic susceptibility of a substance, one first uses a standard compound whose χ is known, for calibration of a particular Gouy apparatus. This essentially involves determining the change in weight of a given mass of standard in the presence and absence of the magnetic field. With a knowledge of the change in weight produced by a given magnetic susceptibility, χ , a measurement of the change in weight, in and out of the field, of an unknown compound will allow its χ to be calculated.

EXPERIMENTAL PROCEDURE

Scratch a horizontal line on the Gouy tube about 2 cm from the top, if this has not already been done. For the weight measurements, the tube should always be filled to this line. To remove paramagnetic impurities, clean the tube with a cleaning solution (not $K_2Cr_2O_7-H_2SO_4$; see Introduction, p. 6), thoroughly rinse with water and acetone, and dry in the oven. Do not wipe the tube with a dry towel; this gives the tube a static charge that significantly affects the weighings.

Weigh the empty tube on the chain with the field off. With the field on, weigh the tube again. Although "pure glass" is diamagnetic, paramagnetic impurities may cause the tube to be attracted by the field rather than repelled. The difference between these two weighings, "on" - "off", is δ and will be used to correct for the magnetism of the tube when the sample is weighed in the tube. It is necessary to maintain the same magnetic field for all measurements when the field is on. With electromagnets this will require a constant electric current (and therefore field). Since the current will decrease as the coils begin to heat, the current will have to be adjusted frequently. Current regulators are available that will conveniently provide this control.

To correct for the magnetism of air when it is displaced by the sample, the volume occupied by the air must be determined. This is done by weighing the tube filled to the line with water. From the weight and known density, d , of water at the existing temperature, the volume (V) may be calculated. The volume susceptibility of air is 0.029×10^{-6} per ml.

To determine the calibration constant for the apparatus, the dry tube is filled to the line with the solid standard. The largest source of

error in the Gouy method is inhomogeneously packed sample tubes. To minimize this problem, the sample should be finely powdered (use a mortar and pestle) and introduced into the tube in small portions. After each addition, firmly tap the tube on a hard surface. Careful packing of the tube will require 20 to 30 minutes. Weigh the tube with the magnet off and again with the magnet on, using the same current as used previously. After each weighing of a solid sample with the field on, measure the temperature between the poles of the magnet. The difference in these two weighings ("on" - "off") is designated Δ and is a measure of the magnetic susceptibility of both the sample and the tube. Knowing the magnetic susceptibility per gram, χ , of the standard, the mass (in grams) of the standard (m), and the values of δ , Δ , and V , the calibration constant, β , may be calculated:

$$\chi(m) = (0.029 \times 10^{-6})V - \beta(\Delta - \delta) \quad (10)$$

For calibration, either $HgCo(NCS)_4$ or $[Ni(en)_3]S_2O_8$, where $en = NiH_2C_2H_4CH_2NH_2$, has proved to be very satisfactory. They may be prepared easily in high purity, are stable, are not hygroscopic, and pack very well. The susceptibility per gram, χ , of $HgCo(NCS)_4$ at 20° is $\chi_{300} = 16.44 \times 10^{-6}$ cgs (± 0.5 per cent). The susceptibility obeys the Curie-Weiss Law with $\theta = 10^\circ$. The relatively high susceptibility of this compound sometimes causes the sample tube to cling to one of the poles of the magnet. This can occasionally be avoided by carefully positioning the sample tube midway between the poles. On the other hand, $[Ni(en)_3]S_2O_8$ is rarely drawn toward a pole because of its lower susceptibility. Its value of χ_{300} is 11.03×10^{-6} cgs (± 1 per cent). It, too, obeys the Curie-Weiss Law with $\theta = -43^\circ$.

After the value of β is obtained, the sample tube is cleaned and dried. The same procedure is repeated for the determination of the unknown, $[Cr(NH_3)_6](NO_3)_3$. The empty tube is weighed with the field on and off to obtain δ . Then the tube is carefully packed with finely powdered $[Cr(NH_3)_6](NO_3)_3$. The filled tube is then weighed with the field on and off to obtain Δ . These measurements will then permit the calculation of χ and the molar susceptibility, χ_M .

Summarized below are the weighings that must be made first on the standard and then on the unknown:

- A. Weight of empty tube, field off ————g
- B. Weight of empty tube, field on ————g
- C. Weight of tube filled to line with water, field off ————g
- D. Weight of tube filled to line with solid, field off ————g
- E. Weight of tube filled to line with solid, field on (measure temperature) ————g

As given, these weights are related to the terms in equation (10) in the following manner:

$$V = (C - A)/d, \text{ where } d \text{ is density of water (g/ml) at ambient temperature}$$

$$\delta = B - A$$

$$\Delta = F - D$$

$$m = D - A$$

To determine the reproducibility of your value, the evaluation of χ for $[\text{Cr}(\text{NH}_3)_6](\text{NO}_3)_3$ should be repeated at least one more time by emptying and repacking the tube and then making the necessary weightings.

REPORT

Include the following:

1. The value of χ for $\text{HgCo}(\text{NCS})_4$ and/or $[\text{Ni}(\text{en})_3]\text{S}_2\text{O}_8$ at the temperature of your weightings.
2. Data necessary to calculate β , χ , χ_M , χ_{SP} and μ_{eff} for $[\text{Cr}(\text{NH}_3)_6](\text{NO}_3)_3$.
3. Sample calculation of each of the above quantities, as well as the number of unpaired electrons in the complex.
4. Account for your experimentally determined number of unpaired electrons in terms of the electronic structure of the complex.

QUESTIONS

1. Account for the increase in the values of the diamagnetic corrections (Table 5-1) in going from the top to the bottom of any group in the periodic table.
2. Using the known values of χ for $\text{HgCo}(\text{NCS})_4$ and $[\text{Ni}(\text{en})_3]\text{S}_2\text{O}_8$, attempt the calculation of the number of unpaired electrons in each of these complexes, and interpret the results in terms of the electronic structures of the complexes.
3. An $\text{Fe}(\text{II})$ complex of the formula FeL_6^{2+} might have either square pyramidal or trigonal bipyramidal geometry. Would it be possible to differentiate between these geometries on the basis of a magnetic susceptibility measurement of the complex? How would the number of unpaired electrons differ in the two geometries? (Use crystal field theory to show splitting of d orbitals for the two structures.)
4. Rationalize equation (10).
5. The volume susceptibility of air is listed as 0.029×10^{-6} per ml. What causes its susceptibility to be positive?
6. If the standard were packed more compactly than the unknown, but both were packed evenly in the tube, how would this affect the calculated value of χ ?

7. Why should $\text{K}_2\text{C}_2\text{O}_7\text{--H}_2\text{SO}_4$ cleaning solution not be used to clean Gouy tubes?

8. A student measured the effective magnetic moment of $\text{NiCl}_2 \cdot 6\text{H}_2\text{O}$ and obtained a value of 3.2 B.M. A value of 0.0 B.M. was obtained for $\text{Ni}(\text{CN})_5^{3-}$. Rationalize these values in terms of the structures and crystal field splitting diagrams of these ions.

9. When the bottom of the tube is at the center of the pole faces, a Gouy tube containing $\text{HgCo}(\text{NCS})_4$ has a certain weight in the magnetic field. If the magnet is raised so that the bottom of the tube is slightly below the center of the pole faces, will the observed weight be more or less than in the original measurement? Explain.*

INDEPENDENT STUDIES

A. Determine the effective magnetic moment and number of unpaired electrons in $[\text{Co}(\text{NH}_3)_2\text{Cl}][\text{Cl}_2 \cdot \text{CrCl}_3 \cdot \text{K}_3[\text{Fe}(\text{C}_2\text{O}_4)_3] \cdot 3 \text{H}_2\text{O}$ (R. C. Johnson, *J. Chem. Educ.*, **47**, 702 (1970)), $\text{K}_3[\text{Fe}(\text{CN})_6]$, or other compounds.

B. Using the Evans solution nmr method, determine the effective magnetic moment and number of unpaired electrons in $\text{K}_3[\text{Fe}(\text{CN})_6]$ and other iron complexes (J. L. Deutsch and S. M. Polling, *J. Chem. Educ.*, **46**, 167 (1969)); J. Lohger and R. Scheffold, *J. Chem. Educ.*, **49**, 646 (1972); D. Ostfeld and I. A. Cohen, *J. Chem. Educ.*, **49**, 829 (1972).)

C. Determine the electron spin resonance spectrum of $\text{K}_3[\text{Mo}(\text{CN})_6]$ or $\text{K}_3[\text{W}(\text{CN})_6]$ (B. R. McGarvey, *Inorg. Chem.*, **5**, 476 (1966).) For general discussion of the electron spin resonance technique, see the following references: A. Abragam and B. Bleaney, *Electron Paramagnetic Resonance of Transition Ions*, Oxford University Press, London, 1970; R. S. Drago, *Physical Methods in Inorganic Chemistry*, Reinhold Publishing Corp., New York, 1965, Chapter 10; B. A. Goodman and J. B. Raynor, *Adv. Inorg. Chem. Radiochem.*, **13**, 135 (1970); B. R. McGarvey, *Transition Metal Chemistry*, R. I. Carlin, ed., Marcel Dekker, New York, 1966, p. 89.

REFERENCES

Magnetochemistry—Techniques and Interpretation

- U. G. Duncanson, *J. Chem. Educ.*, **44**, 142 (1967).
 A. Farnshaw, *Introduction to Magnetochemistry*, Academic Press, London, 1968.
 B. N. Figgis and J. Lewis in *Technique of Inorganic Chemistry*, Vol. 11, H. R. Jonassen and A. Weissberger, eds., Interscience Publishers, New York, 1965, p. 137.
 B. N. Figgis and J. Lewis in *Modern Coordination Chemistry*, J. Lewis and R. G. Wilkins, eds., Interscience Publishers, New York, 1960, p. 460.
 P. W. Selwood, *Magnetochemistry*, Interscience Publishers, New York, 1956.

The Study of a Cobalt Complex—A Laboratory Project

James H. Loehlin, Stephen B. Kahl, and Jeanne A. Darlington,
Wellesley College, Wellesley, MA 02181

Over the past 10 to 15 years at Wellesley, various members of the department have developed a laboratory program which includes an eight-week project that has been very successful. The original experiments which served as the basis for this program were reported by Sebera (1) and Wilson (2). The project is undertaken by all of the students in our second semester introductory course, ordinarily utilizing the first five and last three of the twelve laboratory sessions, although other schedules have been used from time to time.

The laboratory involves synthesis followed by qualitative and quantitative analysis of either $\text{Co}(\text{NH}_3)_5\text{Cl}_3$ or $\text{Co}(\text{NH}_3)_6\text{Cl}_3$. The students are told that the complexes contain cobalt, ammonia, and chlorine but are given no other information about them. During the middle of the semester the students work for four weeks on a series of experiments designed to present standard qualitative analysis techniques and procedures which are not directly related to the cobalt project. Later, they perform three project experiments which are designed to give additional insights into the structure and properties of the complexes prepared earlier. Although any part of this project can stand by itself as a worthwhile experiment, the value of the entire project is much greater than the sum of its parts. The several experiments complement each other very well.

The preparation, purification, and analysis give the student experience with synthesis, recrystallization, gravimetric analysis, and volumetric analysis using both acid-base and redox methods. In addition, she gains experience with qualitative analysis, determining an empirical formula from percent composition by weight and obtaining theoretical and actual yields. At the end of the semester, the student determines the number of unpaired electrons in her compound, the effect of ligands on splitting of the *d* subshell and the actual charge on the complex cation. These final three experiments, which we have developed, are explained in detail below.

At the beginning of each laboratory period, a short discussion is presented, dealing with such topics as methods of purification and crystallization, laboratory instruments and their use, distillation techniques, and volumetric and gravimetric analyses in addition to comments about specific details of the day's experiment. Before each of the analyses, students are assigned a calculation based on hypothetical student data

from the appropriate analysis of $(\text{Co}(\text{NH}_3)_4\text{Cl}_2)\text{Cl}$. Since each student uses the same data for this, it is easy for the instructor to spot errors in her calculation or understanding.

One other benefit accrues during the first five-week period. In lecture, the first several weeks are given to discussions of atomic structure, periodic properties, and theories of the chemical bond. It is most difficult to find suitable introductory level laboratory exercises to illustrate these topics. Thus, the quantitative and qualitative aspects of this complex ion project serve to reinforce and clarify the techniques introduced during the previous semester's laboratory work and to provide a link between the two semesters.

The laboratory work in the semester is concluded with a series of three experiments giving specific structural information about the complex which the student has prepared along with that of other complexes. Most of our students report that these particular experiments are very useful in bringing their knowledge and understanding together. In addition, they seem to enjoy them.

Weekly Outline of the Cobalt Complex Project¹

Week 1 Preparation of a Complex

The cobalt chloropentaammine and hexaammine chlorides are prepared from $\text{Co}(\text{H}_2\text{O})_6\text{Cl}_2$ by the methods of Hynes, Yanowsky, and Shiller (3), and Palmer (4), respectively. These procedures have been published in THIS JOURNAL with minor modifications by Sebera (1) and Wilson (2). The crude product is saved for purification the following week.

Week 2 Purification and Qualitative Analysis

Sebera's and Palmer's procedures are used for the purification and recrystallation of the samples which are then oven dried at 110°C before weighing. A small sample of the wet product is usually removed for use in the qualitative analysis. The presence of NH_3 is shown by the color change of moist litmus paper above a basic solution of the complex. Chloride's presence is shown by AgCl precipitation. Cobalt is first shown to be present by the formation of a red-brown 2-ni-

¹ A packet of information containing detailed instructions used by our students, lab prep lists, and a listing of "Magnetic," a BASIC program to assist students in obtaining μ from experimental weights, may be obtained by writing JHL.

troso-1-naphthol complex. That it is no longer Co^{2+} is demonstrated by a negative thiocyanate complex test which is compared with the same test using a solution of the Co^{2+} starting material.

Week 3 Ammonia Analysis


The ammonia analysis is carried out by driving the NH_3 gas from the complex by heating it in a standard distillation apparatus following the addition of 10% NaOH. The NH_3 is bubbled through a known volume of standardized HCl which is then diluted in a volumetric flask, and aliquots are titrated to determine the amount remaining. The difference between the original amount of acid and that remaining is used to calculate the percent by weight of ammonia in the original sample.

Week 4 Cobalt Analysis

The residue of Co_2O_3 remaining in the distilling flask from the previous week's experiment is used for the Co analysis. An excess of KI/HCl is added and the solution containing the dark brown I_3^- ion is transferred to a volumetric flask, diluted to volume, and aliquots titrated with standardized thiosulfate solution. The titration is carried out to a light golden color before fresh starch indicator is used to determine the final endpoint. Very careful student titrations can turn blue again within a few minutes due to air oxidation of the I^- in the absence of excess thiosulfate. This serves as a nice check on the extent of over-titration.

Week 5 Chloride Analysis

Gravimetric determination of chloride using AgCl is performed in duplicate by each student. The samples are prepared, and the AgNO_3 is added during week 4. They are then left to digest for a week before filtering, drying to constant mass, and performing the calculations. The results of all her analyses are then used to determine the empirical

Detail of Gouy balance showing a sample both in and out of the magnetic field.

formula and the % yield. We have chosen to continue using this method in spite of the drastic increase in the price of silver. We have, however, begun to be very careful to save excess reagents and the students' precipitates for recovery of the silver.

Typical analytical results of groups of students show standard deviations in the cobalt and ammonia analyses of about 1.5%. Since all of our students have previously been exposed to volumetric analysis in the first semester of the course, we generally require those who do not obtain 1% precision in the ammonia and cobalt determinations to titrate additional aliquots. Often the results of the ammonia analyses are slightly low, probably from the loss of ammonia in the distillation, while those of the cobalt analyses are slightly high as a result of the slow appearance of the starch endpoint. The chloride analyses are usually somewhat better than the others, showing standard deviations of around 1%.

Week 6 Summary Report

At the end of this five-week period, each student is required to submit a report summarizing her quantitative and qualitative results. This report includes calculation of an empirical formula for the complex, percentage yield, and a discussion of possible sources of error. When the data do not support a clear choice of empirical formula, a student is encouraged in her sources-of-error discussion to suggest the most plausible choice. We suggest that she consider the implications of such factors as the difference between the sum of her % composition values and 100% and the charge balance among the components of the compound. Throughout the entire process, we try to emphasize the parallels between her experiments and those of a research chemist faced with a similar problem.

Final Weeks

After the above five weeks spent on the cobalt complex project, our students spend four weeks working on various qualitative analysis experiments. During this period transition metal chemistry is introduced in lecture. Rather than choosing some new and unfamiliar systems for illustrating transition metal properties, the last three weeks of laboratory are spent in further study of the cobalt complexes. Students perform the following three experiments in no particular order: Determination of the Number of Unpaired Electrons by Magnetic Susceptibility, Determination of Charge on the Complex Ion by Ion Exchange Chromatography, and Visible Absorption Spectra and the Spectrochemical Series.

Magnetic Susceptibility—Determination of Number of Unpaired Electrons

Introduction and Theory

We have designed a magnetic susceptibility experiment using the Gouy method which allows our students to calculate the number of unpaired d -electrons in a series of transition metal complexes. This is then correlated with the electron configuration of the metal atom. Although several reports concerning the use of the Gouy method for magnetic susceptibility determination have appeared in THIS JOURNAL (5-11), we believe this is the first report of its quantitative use in an introductory course and is the most complete in terms of the spectrum of first row transition metal complexes used. With the possible exception of the magnet, the apparatus used is relatively inexpensive or available in most laboratories.

Samples of complexes with equal transition metal concentration and identical geometry are weighed in and out of the magnetic field gradient and the weight differences, ΔW_{meas} , are calculated ($\Delta W_{\text{meas}} = W_{\text{in}} - W_{\text{out}}$). All paramagnetic samples are corrected for diamagnetic effects using the average weight difference, $\overline{\Delta W}_{\text{dia}}$, of five diamagnetic samples which contain no transition metal atoms. ($\overline{\Delta W}_{\text{dia}} = \frac{1}{5}(\text{Sum}(\Delta W_{\text{meas}}))$). Note that $\overline{\Delta W}_{\text{dia}}$ is negative! The difference in weight due solely to paramagnetic effects, ΔW_{para} , is then

$$\Delta W_{\text{para}} = \Delta W_{\text{meas}} - \overline{\Delta W}_{\text{dia}} \quad (1)$$

This difference in weight is related to the total magnetic moment, μ , by the equation

$$\mu^2 = k \Delta W_{\text{para}} \quad (2)$$

where k is the proportionality constant for the apparatus. The

value of k may be determined from measurements on a sample whose susceptibility is known. The susceptibility of a transition metal ion arises from two factors,

$$\mu = \mu_s + \mu_m \quad (3)$$

where, μ_s is the spin-only magnetic moment, and μ_m is the magnetic moment due to spin-orbit coupling. The spin-only magnetic moment, μ_s , may be calculated by

$$\mu_s = \sqrt{n(n+2)} \text{ Bohr magnetons} \quad (4)$$

where n = number of unpaired electrons.

If care is taken to keep metal concentrations, magnetic field and sample geometry relative to the field constant, k will be the same for all samples. For high-spin d^5 octahedral complexes, the environment should differ little from spherical symmetry where μ_m is equal to zero. Assuming a zero value for μ_m , we can calculate μ from eqns. (4) and (3) and then use this value along with the experimentally measured ΔW_{para} in eqn. (2) to evaluate k . Once k has been determined, eqn. (2) may be used to calculate μ for all other paramagnetic samples. By comparing experimental values of μ with all the possible spin-only moments calculated from eqn. (4), the number of unpaired electrons in any sample may be found by inspection, since $\mu \geq \mu_s$.

Experimental

Our apparatus consists of an electromagnet with adjustable gap and 5 in. poles tapered to 4 in. at the face (Alpha Scientific Labs, Model 500) with its associated power supply (AL 7500) and regulator (7500 R). These models are no longer manufactured, but any magnet producing a field of several kilogauss in a gap about 2 or 3 in. in diameter and 1 in. wide should be satisfactory. The magnet is mounted on the lower shelf of a table (see the figure). A Mettler H4 balance is mounted on the table directly above the magnet gap with a chain attached to the pan passing through the base of the balance and the table. Samples are hung at the end of the chain for in-field weighings and from the balance pan hook for out-of-field weighings. The entire apparatus must be carefully levelled to ensure that, in the absence of a field, the sample hangs in the center of the gap and is parallel with its sides. A guide is used to keep the samples from swinging against the pole faces when the magnet is on. The bottom of the sample tube hangs in the center of the magnet gap where the field strength is at its maximum. We have also enclosed the chain and magnet gap in a rigid plastic enclosure to minimize the effect of drafts on the weighings.

Students are provided with 15 ml of 1.00 M aqueous solutions of the soluble compounds sealed in matched 18 mm test tubes. (Three of the complexes, the two synthesized by the students and $K_4Fe(CN)_6$ are prepared using a corresponding amount of complex but are not completely soluble at room temperature so consist of a saturated solution with crystals in the bottom. The average concentration of each 15 ml sample is still 1.00 M. The resulting inhomogeneity might lead to inaccuracies in determination of the number of unpaired electrons except that these three are all diamagnetic. The chloropentaammine is hydrolyzed slowly to the aquopentaammine complex in aqueous solution, but this does not interfere with this experiment.) Copper wire hangers are attached to the tubes with epoxy cement so that the sample geometry relative to the magnet gap is identical in all samples. A current of seven amps gives a field strength of ~ 8000 gauss at the gap center. This field produces a weight difference of approximately 1 g for the Mn^{2+} or Fe^{3+} samples.

Each student in a laboratory section weighs two samples in and out of the field. She is later given the complete set of data for the section. At least one set of measurements on each sample (the twelve solutions shown in Table 1 and 1.0 M solutions of $Ca(NO_3)_2$, $NaCl$, NH_4Cl , KBr , and H_2SO_4 for determination of the diamagnetic correction) are given. Following the calculation of μ either directly or with the aid of the computer program, "Magnetic," she interprets the result in terms of the number of unpaired electrons for each transition metal atom and compares this with the number of d electrons and the splitting of the d orbitals.

Results and Discussion

The results from this experiment are shown in Table 1. The derivation of n , the number of unpaired electrons, is unambiguous using the simplified theory presented here but does

require some care on the part of students in interpreting their data. The value for Cr(III) is slightly below the spin-only value of 3.87 Bohr magnetons, which is within the limits we suggest to students for uncertainties resulting from theoretical approximations and experimental error. The value for Co(II), on the other hand, is larger than the spin-only value for four unpaired electrons, 4.90 Bohr magnetons, but since there is an odd number of electrons in the ion, n must also be odd and hence equal to three.

Several other points should be noted. First, agreement between observed and reported values is excellent. The resulting values are reproducible within about 3% from one laboratory period to another and for different samples of the same solution. Second, nine of the fourteen possible configurations of ions with d^1 through d^{10} configurations are represented. Octahedral high- and low-spin d^4 and low-spin d^7 complexes of simple ligands are generally unstable in solution, while samples of d^1 and d^2 ions, e.g., $(Ti(H_2O)_6)^{3+}$ and $(V(H_2O)_6)^{3+}$, respectively, could be included, but would require using inert atmosphere conditions. We have chosen to omit them from this experiment, since all of the other samples are either commercially available or can be synthesized easily. Third, for d^5 and d^6 configurations, the effect of ligand field can be seen clearly.

Absorption Spectra—Determination of Spectrochemical Series

Introduction and Theory

This experiment is designed to make a direct determination of the spectrochemical series for several ligands. Students measure the visible absorption spectra of seven Co(III) complexes. These spectra are typical of octahedral low-spin d^6 complexes, which show transitions from the ${}^1A_{1g}$ ground state to ${}^1T_{1g}$ and ${}^1T_{2g}$ excited states. The $A_{1g} \rightarrow T_{1g}$ transition is lower in energy, and represents the magnitude of the crystal field splitting energy, $10Dq$ or Δ_0 . Since the second transition, $A_{1g} \rightarrow T_{2g}$, has a maximum absorbance well beyond the high energy end of the 450–600 nm region examined by the students, only the lower energy side of this peak is observed. Only

Table 1. Magnetic Susceptibilities of First Row Transition Metals

| Metal ion | Solution prep. from | μ_{exp} Bohr Magnetons | μ_{th} | [From Figgis and Lewis (12) Reference Compound | page |
|-----------------|---|---|-------------------|--|------|
| d^3 Cr(III) | Cr(NO ₃) ₃ | 3.85 | 3.83 | CrCl ₃ ·6H ₂ O | 128 |
| d^5 Mn(II) | MnSO ₄ ·H ₂ O | 5.91 ^a | 5.86 | MnSO ₄ ·4H ₂ O | 155 |
| | FeCl ₃ ·6H ₂ O | 5.92 ^a | 5.95 | FeCl ₃ ·6H ₂ O | 169 |
| | K ₃ Fe(CN) ₆ | 2.33 | 2.25 | K ₃ Fe(CN) ₆ | 172 |
| d^6 Fe(II) | FeSO ₄ ·7H ₂ O | 5.46 | 5.47 | (NH ₄) ₂ Fe(SO ₄) ₂ ·6H ₂ O | 177 |
| | K ₄ Fe(CN) ₆ ·3H ₂ O | 0 | — | | |
| | Co(III) | Co(NH ₃) ₆ Cl ₃ | 0 | — | |
| | Co(NH ₃) ₆ Cl ₃ | 0 | — | | |
| d^7 Co(II) | Co(NO ₃) ₂ ·3H ₂ O | 4.92 | 4.84 | Co(NO ₃) ₂ ·6H ₂ O | 188 |
| d^8 Ni(II) | Ni(NO ₃) ₂ ·6H ₂ O | 3.24 | 3.38 | Ni(H ₂ O) ₆ Cl ₂ | 200 |
| d^9 Cu(II) | CuCl ₂ | 1.99 | 1.95 | CuSO ₄ ·5H ₂ O | 213 |
| d^{10} Zn(II) | Zn(NO ₃) ₂ ·6H ₂ O | 0 | — | | |

^a Value assumed to be 5.916 in calculating k .

Table 2. Absorbance Maxima for Cobalt Complexes

| Complex Ion | Differentiating Ligand | λ_{max} (nm) | Prepared Compound |
|--|-------------------------------|-----------------------------|---|
| (Co(en) ₃) ³⁺ | en | 467 | Co(en) ₃ Cl ₃ |
| (Co(NH ₃) ₆) ³⁺ | NH ₃ | 475 | Co(NH ₃) ₆ Cl ₃ |
| (Co(NH ₃) ₅ H ₂ O) ³⁺ | H ₂ O | 495 | Co(NH ₃) ₅ H ₂ O(NO ₃) ₃ |
| (Co(NH ₃) ₅ OCO ₂) ⁺ | CO ₃ ²⁻ | 510 | Co(NH ₃) ₅ (CO ₃)NO ₃ |
| (Co(NH ₃) ₅ F) ²⁺ | F ⁻ | 515 | Co(NH ₃) ₅ F(NO ₃) ₂ |
| (Co(NH ₃) ₅ Cl) ²⁺ | Cl ⁻ | 534 | Co(NH ₃) ₅ Cl ₃ |
| (Co(NH ₃) ₅ Br) ²⁺ | Br ⁻ | 552 | Co(NH ₃) ₅ Br ₃ |

Results and Discussion

Precision and accuracy of the experimental measurements are excellent, generally falling within 1–2% of each other and of the true value. The results of the experiment show the +2 or +3 complex ion charge, depending on the complex used. Since the charge on the cobalt metal ion itself has been found previously by qualitative analysis, students can infer that one of

the hexaammine cobalt(III) complex has true octahedral symmetry. The symmetry reduction for the tris(en) complex and the ligand-pentaammine mixed ligand complexes does not result in observable splitting of the T_{1g} excited state but does broaden slightly the transition peak area.

Experimental

the hexaammine cobalt(III) complex has true octahedral symmetry. The symmetry reduction for the tris(en) complex and the ligand-pentaammine mixed ligand complexes does not result in observable splitting of the T_{1g} excited state but does broaden slightly the transition peak area.

Experimental

Spectra are measured using Bausch and Lomb Spectronic 20 s at 10 nm intervals over the 450–600 nm range. To get a more accurate measure of λ_{\max} , the absorbances are then remeasured at 5 nm intervals in the peak region. The single beam design of these instruments makes it necessary to re-standardize with a water blank before each reading. Pairs of students, each of whom initially prepared a different complex, work together. In addition to their two complexes, a sample of one of the remaining five complexes shown in Table 2 is supplied to each pair. These compounds were synthesized using procedures given by Schlesinger (13). Each pair is then asked to prepare 10 ml of a 0.010 molar solution of each compound (5 millimolar for the bromo and carbonato complexes) given the formula weight. They report the absorbances of the three samples at 500 nm to the instructor as a check on solution preparation before proceeding to measure and graph the entire set of spectra. The experimentally determined values of λ_{\max} are shared around the class. A full set of seven λ_{\max} 's is thus available to each student.

Results and Discussion

Our experimental data, given in Table 2, show that the peak maxima shift in agreement with reported spectrochemical series for the differentiating ligands in these complexes. Students are encouraged to use the Planck equation to convert the wavelength values to energy units to give a direct evaluation of the splitting energy. The spectrochemical series obtained is compared with the one in the students' text which usually is missing some of the ligands studied. They, thus, not only verify the text information and see how it was obtained but also expand it to include new ligands.

Ion Exchange—Determination of Charge on a Complex Cation

Introduction and Theory

The charge on the complex ion is arrived at by percolating an aqueous solution of a known weight of compound through a bed of acidified cation exchange resin followed by titration of the acidic eluent. To maintain electroneutrality, each $n+$ charged ion absorbed by the resin releases n H^+ ions which are collected in the eluent and titrated. Since the empirical formula of the unknown complex has been determined previously, the charge on the complex ion may be calculated easily from the moles of acid and of compound. This experiment is an adaptation of one given by Bell (14) where similar ion-exchange methods are used to determine the amount of cobalt in $Co(en)_2Cl_3$.

Experimental

We have used Dowex-20 resin whose structure consists of a styrene-divinylbenzene copolymer with sulfonate groups as the active exchange sites. The protonated form of the resin is generated by soaking in 1.5 M HCl for a few minutes. It is then rinsed free of excess H^+ with distilled water before the column is packed. A slurry of resin (12 g) in water is poured into a buret previously fitted with a small plug of glass wool. (We have employed old pinch-clamp burets which were no longer being used here.) The resin settles to form a column approximately 6 in. deep. This is sufficient to allow the column to be used for at least 3 or 4 samples, although most students get suitable data from duplicate runs. The packed column is rinsed with distilled water until the eluent pH is the same as the water. To prevent air entrapment in the column, the liquid level must never be allowed to drop below that of the top of the resin bed.

A weighed sample of 100–150 mg of complex is dissolved in 50 ml of water and percolated through the column with distilled water. The eluent is tested with pH paper until its pH matches that of distilled water to ensure collection of all protons generated. The amount of H^+ released is determined by titration with standardized 0.05 M NaOH to a phenolphthalein endpoint. The procedure is repeated with a second sample of complex on the same column. The resin may be collected and regenerated for use by later classes by soaking overnight in 6 M HCl and rinsing with distilled water.

Results and Discussion

Precision and accuracy of the student determinations are excellent, generally falling within 1–2% of each other and of the true value. The results of the experiment show the +2 or +3 complex ion charge, depending on the complex used. Since the charge on the cobalt metal ion itself has been found previously by qualitative analysis, students can infer that one of the three chloride ions in $Co(NH_3)_6Cl_3$ is part of the octahedral complex ion, or that only neutral ammonia molecules are bound to the cobalt in $Co(NH_3)_6Cl_3$.

Summary and Conclusions

These three final experiments tie together the lecture material on the structure of transition metal complexes with first-hand laboratory experience in obtaining relatively sophisticated structural information. With the empirical formula of her compound in hand and considerable information about octahedral transition metal complexes from class, the student can use the experimentally determined complex ion charge to infer the nature of the ligands. She is then asked to describe the nature of the bonding between the transition metal and ligands in comparison with that of the nonligand chlorides with the complex ion.

The ligands of her compound are next placed into the spectrochemical series which serves to relate the d electron energy level discussion from class with experimentally determined transitions. We have found it useful also to discuss one of the experimental spectra and how it relates to the observed color of the complex. The effect of the cyanide ligand on the magnetic susceptibility of the hexacyanoferrate(II) and (III) complexes allows correlation of one more ligand with the spectrochemical series.

The direct determination of the number of unpaired electrons using the magnetic susceptibility experiment brings together the concepts of electronic structure, Aufbau order, and the perturbing effects of ligands on a degenerate energy level. Asking the students to predict the color of a zinc ion solution during the final oral discussion ties many of these ideas together. The magnetic susceptibility sample may then be brought out to compare with the predictions.

The cobalt complex project provides the student with a wide variety of experiments, each of which adds to her knowledge of the compound synthesized in the first laboratory period. At the same time, each week's experiment uses a different experimental approach, thereby broadening the student's range of laboratory procedures.

Students find these experiments interesting and challenging, and almost universally they agree that the project orientation serves to unify their laboratory experience throughout the semester. We have also found that the last three experiments herein described bring life and vitality to the lecture and text discussions of transition metal chemistry and bonding.

Literature Cited

- (1) Sebera, D. K., *J. CHEM. EDUC.*, **40**, 476 (1963).
- (2) Wilson, L. R., *J. CHEM. EDUC.*, **46**, 447 (1969).
- (3) Hynes, W. A., Yarowski, L. K., and Shiller, M., *J. Amer. Chem. Soc.*, **60**, 3053 (1938).
- (4) Palmer, W. G., "Experimental Inorganic Chemistry," Cambridge Univ. Press, Cambridge, England, 1954, pp. 530, 534.
- (5) Eaton, S. S. and Eaton, G. R., *J. CHEM. EDUC.*, **56**, 170 (1979).
- (6) Viswanadham, P., *J. CHEM. EDUC.*, **55**, 54 (1978).
- (7) Bobonich, H. M., *J. CHEM. EDUC.*, **42**, 495 (1965).
- (8) Brubacher, L. J. and Stafford, F. E., *J. CHEM. EDUC.*, **39**, 574 (1962).
- (9) Kirschner, S., Albinak, M. J. and Bergman, J. G., *J. CHEM. EDUC.*, **39**, 575 (1962).
- (10) Ihrig, J. L. and Calkwell, R. G., *J. CHEM. EDUC.*, **32**, 320 (1955).
- (11) Schuler, R. H., *J. CHEM. EDUC.*, **27**, 591 (1950).
- (12) Figgis, B. N. and Lewis, J., "The Magnetic Properties of Transition Metal Complexes," in "Progress in Inorganic Chemistry," F. A. Cotton, (Editor), Interscience Publishers, J. Wiley & Sons, New York, 1964, Vol. 6, pp. 37–239.
- (13) Schlesinger, G. G., "Inorganic Laboratory Preparations," Chemical Publishing Co., New York, 1962.
- (14) Bell, J. A., "Chemical Principles in Practice," J. A. Bell, (Editor), Addison-Wesley, Reading, MA, 1967, pp. 47–50.

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NMR Spectroscopy Using a Chiral Lanthanide Shift Reagent To Assess the Optical Purity of 1-Phenylethylamine

Tito Viswanathan and Alan Toland

University of Arkansas at Little Rock, Little Rock, AR 72204

To obtain the optical purity of a sample using polarimetry requires the value of the absolute rotation of the pure enantiomer. However, this is a condition that cannot always be met. The separation of a racemic mixture into its constituent enantiomers (i.e., resolution or optical resolution) and the determination of the optical purity of the resolved material has therefore provided a challenge for chemists. The problem is acute because the only differences between enantiomeric pairs can be subtle.

In contrast, diastereomers (i.e., stereoisomers that are not enantiomers) are distinguishable even in achiral surroundings and exhibit different physical properties, for example, melting point, solubility, and IR and NMR spectra. It thus seems logical that most methods of resolution involve the conversion of enantiomers into diastereomers (1).

The resolution of a racemic form involves reaction with a chiral substance (preferably optically pure), called a resolving agent, to produce a mixture of diastereomers. Although the constituents of a racemic mixture have the same solubility in achiral solvents, the diastereomers produced usually do not. They can usually be separated by crystallization from an achiral solvent.

An experiment to exemplify the use of NMR in determining the extent of resolution is being routinely carried out in our organic laboratory. An enantiomeric mixture of 1-phenylethylamine (α -phenylethylamine, α -methylbenzylamine), is reacted with the chiral reagent (*R,R*)-(+)-tartaric acid to give diastereomeric salts that differ greatly in their solubility in methanol. The enantiomeric reactants have the same structure but different configurations, being mirror images; they have identical physical and chemical properties, except for the direction of rotation of plane-polarized light and the



Figure 1. 200 MHz ^1H -NMR spectrum of the chiral LSR- $\text{Yb}(\text{TFC})_3$ in CDCl_3 .


Figure 2. ^1H -NMR spectrum of (*S*)-(-)-1-phenylamine in CDCl_3 .

interaction with chiral molecules. The diastereomeric products have the same structure but different configurations; they are not mirror images.


The (*R*)-(+)-amine (*R,R*)-(+)-tartrate is soluble in methanol, whereas its diastereomeric companion (*S*)-(-)-amine (*R,R*)-(+)-tartrate is relatively insoluble in methanol and hence is easily separated. The amine salt is converted to the free base by addition of concentrated NaOH and then extracted with ether. The experimental protocol for accomplishing this may be found in standard laboratory texts (e.g., see ref 2).

NMR Spectroscopy


Enantiomeric forms of 1-phenylethylamine cannot be distinguished by ^1H NMR or ^{13}C NMR because the groups attached to the stereocenter are in an enantiotopic environment. However, the chemical shifts of the proton and carbons attached to the stereocenter can be differentially altered so that they become distinguishable in the ^1H and ^{13}C NMR spectra. This is accomplished by using one of several chiral lanthanide shift reagents (LSR's) (3). An extensive review of chiral LSR's in organic chemistry has been published (4).

LSR's are paramagnetic complexes of certain lanthanides, such as europium and ytterbium, with ligands designed to make them soluble in organic solvents. They are also chemically inert and in some cases improve the solubility of compounds in nonaqueous solvents. When added to NMR samples, they coordinate weakly to polar functional groups, such as amines, esters, ketones, and alcohols, and create a strong local magnetic field that produces large chemical shift changes.

Apart from their basic use in attempting to simplify unresolved spectra, shift reagents may be used to determine enantiomeric purity. Shift reagents with chiral ligands are available in optically active form. If these bind a substrate possessing a stereocenter, two diastereomeric complexes can be formed from its enantiomers, which in principle will exhibit different chemical shifts. Examples of a chiral LSR and a nonchiral LSR (both of which are commercially available) are shown below.

Tris(6,6,7,7,8,8,8-heptafluoro-2,2-dimethyl-3,5-octaedionato)europium(III)

Tris[3-(trifluoromethyl)hydroxymethylene]-(+)-camphorato]ytterbium(III)

The addition of the chiral shift reagent Yb(tfc)₃ to racemic 1-phenylethylamine results in the formation of two diastereomeric forms of an acid-base complex, thereby creating a diastereotopic environment for the groups attached to the stereocenter. This leads to distinct resonance peaks (¹H or ¹³C) of the two enantiomeric forms of 1-phenylethylamine. Yb(tfc)₃ resonances themselves occur at high fields and move to still higher fields upon addition of a coordinating substrate (amine).

Following the resolution of the racemic amine, NMR in the presence of the chiral LSR is studied and compared with the spectrum of the racemic sample with the chiral agent as well as the spectrum of a commercially available optically pure (*S*)-(-)-amine sample with the chiral agent.

Results and Discussion

Figure 1 shows the 200-MHz ¹H NMR spectrum of Yb(tfc)₃ in CDCl₃. The various overlapping resonance peaks appear over a broad range denying spectral clarity. The proton spectra of (*S*)-(-)-phenylethylamine (or the *R* isomer or the optically impure +/- mixture for that matter) exhibit sharp peaks shown in Figure 2. Chemical shifts of hydrogens attached to nitrogen appear at 1.5 ppm in this compound. Upon coordination with a LSR this peak experiences the greatest shift, down to 25 ppm. Figure 3 shows the results obtained when Yb(tfc)₃ is added (0.25 mol to 1 mol of substrate). The spectral region between 7.5 and 17.0 ppm shows two set of peaks as compared to those in

Figure 2. (Amine protons that appear far downfield are not shown for reasons of clarity.)

Upon coordination with the chiral LSR the protons of the (*S*)-(+)- and (*R*)-(-) isomers of the amine experience a diastereomeric environment (unlike their previous enantiomeric environment). This causes differential downfield shifts and results in the spectrum shown in Figure 3.

In order to determine which peaks are due to the *S* isomer and which are due to the *R* isomer, a spectra is run of pure *S* enantiomer with the chiral LSR. This helps assign the peaks due to each of the isomers. The results are displayed in Figure 3, where the peaks due to the *S* and the *R* isomers are individually assigned.

The optical purity of a particular sample (such as the optically resolved sample prepared by a student) may be determined by acquiring the NMR with the chiral LSR (same molar ratio as earlier) and determining the relative peak intensities of the α-C-H protons that are well-separated (see Fig. 3). This ratio indicates the molar ratios of the two isomers. Complete resolution will yield NMR peaks due to only one isomer. (A background lecture on stereoisomerism and NMR is required to refresh students memories.)

Some quizzes after the prelab (or postlab) have included questions, such as the following.

- What would be the result of the experiment if (*S*),(*S*)-(-)-tartaric acid were used in place of (*R*),(*R*)-(+)-tartaric acid?
- How would you expect the IR spectrum of either enantiomer of 1-phenylethylamine to compare with that of a mixture?
- What is the basis for the use of a chiral shift reagent to determine the enantiomeric composition of an optically impure 1-phenylethylamine sample?
- Can you think of a possible use for (*S*)-(-)-phenylethylamine?

Literature Cited

1. Jacobus, J.; Raban, M. *J. Chem. Educ.* **1969**, *46*, 351.
2. Ault, A. *Techniques and Experiments for Organic Chemistry*, 5th ed.; Allyn and Bacon: Newton, MA, 1987; p 339.
3. McCreary, M. D.; Lewis, D. W.; Wernick, D. L.; Whitesides, G. M. *J. Amer. Chem. Soc.* **1974**, *96*(4), 1038.
4. Sullivan, G. R. *Top. Stereochem.* **1976**, *10*, 287.


Figure 3. ¹H-NMR spectrum of (+)-1-phenylethylamine + Yb(TFC)₃ (1 M: 0.25 M) in CDCl₃.


- [6] M. A. Alvarez, M. E. García, V. Riera, M. A. Ruiz, C. Bois, Y. Jeannin, *J. Am. Chem. Soc.* **1993**, *115*, 3786.
- [7] The spectroscopic data for **1** indicates that it is isostructural to its molybdenum analogue, the crystal structure of which shows terminal carbonyl ligands. K. A. Azam, A. J. Deeming, M. S. B. Felix, P. A. Bates, M. B. Hursthouse, *Polyhedron* **1988**, *7*, 1793.
- [8] K. M. Nicholas, *Organometallics* **1982**, *1*, 1713.
- [9] M. A. Alvarez, M. E. García, V. Riera, M. A. Ruiz, C. Bois, Y. Jeannin, *Angew. Chem.* **1993**, *105*, 1232; *Angew. Chem., Int. Ed. Engl.* **1993**, *32*, 1156.
- [10] a) Crystal structure determination of **5**-acetone: crystals from acetone/toluene, monoclinic, space group $P2_1$; $a = 16.387(8)$, $b = 11.489(9)$, $c = 10.806(5)$ Å, $\beta = 91.34(3)^\circ$, $Z = 2$, $\rho_{\text{calc}} = 1.84$ g cm $^{-3}$. 3734 independent reflections, of which 1262 with $I_o > 3\sigma I$, were used for calculations, $R = 0.057$, $R_w = 0.067$; b) crystal structure determination of **6**: crystals from CH $_2$ Cl $_2$ /Et $_2$ O, monoclinic, space group $P2_1/c$, $a = 18.164(6)$, $b = 11.355(6)$, $c = 17.859(2)$ Å, $\beta = 97.51(1)^\circ$, $Z = 4$, $\rho_{\text{calc}} = 1.91$ g cm $^{-3}$. 6299 independent reflections, of which 3139 with $I_o > 3\sigma I$, were used for calculations, $R = 0.051$, $R_w = 0.058$. Further details of the crystal structure investigations may be obtained from the Director of the Cambridge Crystallographic Centre, 12 Union Road, GB-Cambridge CB2 1EZ (UK), on quoting the full journal citation.
- [11] a) R. H. Fong, C.-H. Lin, H. Idmoumaz, W. H. Hersh, *Organometallics* **1993**, *12*, 503; b) K. A. Johnson, W. L. Gladfelter, *ibid.* **1990**, *9*, 2101; c) A. A. Aitchison, L. Farrugia, *ibid.* **1987**, *6*, 819.
- [12] D. Barrat, S. J. Davies, G. P. Elliot, J. A. K. Howard, D. B. Lewis, F. G. A. Stone, *J. Organomet. Chem.* **1988**, *355*, 231.
- [13] a) R. L. Miller, R. Toreki, R. E. LaPointe, P. T. Wolczanski, G. D. van Duyne, D. C. Roe, *J. Am. Chem. Soc.* **1993**, *115*, 5570 and references therein; b) W. A. Herrmann, *Angew. Chem.* **1982**, *94*, 118; *Angew. Chem. Int. Ed. Engl.* **1982**, *21*, 117.

A New Multifunctional Heterobimetallic Asymmetric Catalyst for Michael Additions and Tandem Michael–Aldol Reactions**

Takayoshi Arai, Hiroaki Sasai, Kei-ichi Aoe, Kimio Okamura, Tadamasu Date, and Masakatsu Shibasaki*


Although numerous asymmetric catalysts (chemzymes) have been developed utilizing nearly all the elements in the periodic table, the abilities of chemzymes have, in general, proven inferior to those of enzymes, probably because chemzymes have only one function with which to control asymmetric reactions. It is possible that multifunctional heterobimetallic asymmetric catalysts could compensate for the shortcomings of chemzymes by providing positive cooperation between different metals. Recently, for catalytic asymmetric carbon–carbon bond-forming reactions we introduced a new class of heterobimetallic complexes that consist of a rare earth metal center, an alkali metal center (Li, Na, K), and 2,2'-dihydroxy-1,1'-binaphthyl (BINOL).^[1, 2] For example, the LaLi–BINOL complex is quite effective in catalyzing asymmetric nitroaldol reactions, and LaNa–BINOL is a very efficient catalyst for asymmetric Michael reactions. Mechanistic studies on LaNa–BINOL-catalyzed Michael reactions have revealed that though this complex acts as a base catalyst (ONa), it simultaneously shows Lewis acid character (La), which makes highly enantioselective reac-

tions possible even at room temperature.^[2] This result suggested that extension of the "heterobimetallic concept" to metal centers other than rare earths could lead to heterobimetallic asymmetric catalysts with novel properties. In this paper we report the new heterobimetallic asymmetric catalyst **I** consisting of aluminum, lithium, and BINOL, which makes possible for the first time catalytic asymmetric tandem Michael–aldol reactions.

We began with the development of an amphoteric asymmetric catalyst made with aluminum and an alkali metal.^[3] Initially, in an attempt to prepare the AlLi–BINOL complex **I**, diisobutylaluminum hydride (1.0 M solution in THF) was treated with two equivalents of (*R*)-BINOL in THF at 0 °C for half an hour, and the resulting solution was further treated with one equivalent of BuLi (1.73 M solution in hexane) at 0 °C for half an hour. We were pleased that this new asymmetric catalyst (10 mol%^[4]) was also effective in the Michael reaction of cyclohexenone (**6**) with dibenzyl malonate (**4**), giving **7** in 98% *ee*^[5] and 46% yield at room temperature. In addition, we observed that the same catalyst **I** was prepared more efficiently from LiAlH $_4$ and two equivalents of (*R*)-BINOL (THF, 0 °C, 1 h). This catalyst (10 mol%^[4]) provided high yields of **7** in 99% *ee*, as well as excellent yields of the various Michael adducts **3**, **5**, **9**, and **11** in enantiomeric excesses ranging from 91% to 98% (Table 1).^[2]

Table 1. Catalytic asymmetric Michael additions promoted by Al–M–(*R*)-BINOL complexes (10 mol%).

| No. | Enone | Donor | Product | M | t [h] | Yield [%] | <i>ee</i> [%] |
|-----|----------|-----------|-----------|----|---------|-----------|---------------|
| 1 | 1 | 2 | 3 | Li | 72 | 84 | 91 |
| 2 | 1 | 4 | 5 | Li | 60 | 93 | 91 |
| 3 | 6 | 4 | 7 | Li | 72 | 88 | 99 |
| 4 | 6 | 4 | 7 | Na | 72 | 50 | 98 |
| 5 | 6 | 4 | 7 | K | 72 | 43 | 87 |
| 6 | 6 | 4 | 7 | Ba | 6 | 100 | 84 |
| 7 | 6 | 8 | 9 | Li | 72 | 90 | 93 |
| 8 | 6 | 10 | 11 | Li | 72 | 87 | 95 |

These are the highest enantiomeric excesses of all known catalytic asymmetric Michael reactions.^[2, 6] Although LaLi–BINOL and LaNa–BINOL complement each other in their ability to catalyze asymmetric nitroaldol and Michael reactions,^[1, 2] Al–M–(*R*)-BINOL complexes (M = Li, Na, K, and Ba)^[7] are commonly useful for the catalytic asymmetric Michael reaction (entries 3–6).

With these excellent results in hand, we set out to determine the structure of **I** unequivocally. Although an X-ray quality crystal^[8] of **I** was obtained, it was prone to decomposition during X-ray diffraction studies. However, the complex **II** generated from **I** and cyclohexenone (**6**) in THF provided a crystal that was stable to X-ray crystal analysis.^[9] The structure is shown in

[*] Prof. Dr. M. Shibasaki, Dr. H. Sasai, T. Arai
Faculty of Pharmaceutical Sciences, University of Tokyo
Hongo, Bunkyo-ku, Tokyo 113 (Japan)
Fax: Int. code + (3)5684-5206

K. Aoe, K. Okamura, T. Date
Analytical Chemistry Research Laboratory, Tanabe Seiyaku Co., Ltd., Osaka
(Japan)

[**] This study was financially supported by a Grant-in-Aid for Scientific Research on Priority Areas from the Japanese Ministry of Education, Science, and Culture.


Fig. 1. Crystal structure of $C_{40}H_{24}O_2AlLi/C_6H_8O/(C_4H_8O)_3$ (**II**).

assumed to be **I** or a corresponding dimeric structure.

The small electronegativity value of lithium (1.0) relative to that of aluminum (1.5) suggests that a lithium enolate should be generated from a malonate derivative. However, the role of aluminum in this reaction is not well understood. To study the interaction between the enone and the aluminum atom, an ^{27}Al NMR study was performed.^[12] One broad signal was observed for the $AlLi-(R)$ -BINOL complex **I** at $\delta = 75$ (S_0). Upon addition of three equivalents of cyclohexenone (**6**) to the complex, two additional signals were obtained at $\delta = 40$ (S_1) and 23 (S_2). The S_2 signal became more intense and was a major peak under conditions of the Michael reaction. This is strong evidence for an octahedral arrangement of six magnetically equivalent ligands around the aluminum cation in THF solution, since all other shifts are upfield relative to S_0 .^[13] These NMR studies clearly indicate that the carbonyl groups of the enones are coordinated to the aluminum. Consequently, this $AlLi-(R)$ -BINOL complex **I** probably also acts as a heterobimetallic multifunctional asymmetric catalyst,^[2] making possible efficient Michael reactions even at room temperature.

This aspect of the mechanism suggested that the reaction of a lithium enolate derived from a malonate derivative with an enone would lead to an aluminum enolate intermediate. Further studies were carried out to obtain direct evidence for an aluminum enolate in an effort to learn how to enhance the aluminum enolate's influence on the chemical fate of the Michael adducts. On the basis of the electronegativity of aluminum, the protonation of the aluminum enolate should be slower than that of the corresponding lithium, sodium, or lanthanoid enolates. Could such an Al-enolate then be trapped by an electrophile such as an aldehyde? As was anticipated, the reaction of cyclopentenone (**1**), diethyl methylmalonate (**2**), and 3-phenylpropanal (**12**) in the presence of 10 mol% of the $AlLi-(R)$ -BINOL complex gave the three-component coupling product **13** as a single isomer in 91% *ee* (64% yield; Table 2). Quite interestingly, the use of $LaLi-BINOL$, $LaNa-BINOL$, or a lithium-free $La-BINOL$ complex^[6] gave very unsatisfactory results in


Table 2. Tandem Michael-aldol reactions.

| Catalyst | 3 | | 13 | |
|---------------------|-----------|---------------|-----------|---------------|
| | Yield [%] | <i>ee</i> [%] | Yield [%] | <i>ee</i> [%] |
| $AlLi-(R)$ -BINOL | 7 | 90 | 64 | 91 |
| $LaLi-(R)$ -BINOL | 46 | 3 | 30 | – [a] |
| $LaNa-(R)$ -BINOL | 73 | 86 | trace | – |
| $La-(R)$ -BINOL [b] | 57 | 83 | trace | – |


[a] Nonseparable mixture. [b] See Ref. [6].

Figure 1.^[10] The complex has tetrahedral geometry around the aluminum center with average Al–O distances of 1.75 Å. The long Li–O(1) distance of 1.96 Å is indicative of the ionic character of the interaction between Li^+ and $[Al(binol)_2]^-$.^[11] The structure of the $AlLi-BINOL$ complex is thus

the three-component coupling reaction. This is probably because the lanthanum enolate is more reactive toward the acidic proton.^[14] Moreover, the use of benzaldehyde (**14**) instead of **12** gave the three-component coupling product **15** in 82% yield. Although **15** was a mixture of diastereomers, its oxidation gave the corresponding diketone derivative **16** as a single isomer in 89% *ee*. This is the first example of a catalytic asymmetric tandem Michael-aldol reaction.

We therefore describe the reaction pathway in these three-component coupling reactions as follows (Scheme 1). The reaction of diethyl methylmalonate (**2**) with the $AlLi-(R)$ -BINOL complex gives the corresponding lithium enolate (**III**). Enolate **III** then reacts with cyclopentenone (**1**), which is pre-coordinated to the aluminum, to give the aluminum enolate **IV** enantioselectively. Further reaction of enolate **IV** with aldehyde would lead to an alkoxide (**V**). Although it is unclear whether the

Scheme 1. Possible mechanism for the catalytic asymmetric tandem Michael-aldol reaction promoted by the $AlLi-(R)$ -BINOL complex (**I**).

aluminum or lithium alkoxide is generated, the resulting alkoxide then abstracts a hydrogen atom from an acidic OH group to give the three-component coupling product and regenerate the $AlLi-(R)$ -BINOL complex, thereby completing the catalytic cycle.

In conclusion, we have succeeded in developing highly efficient $Al-M-(R)$ -BINOL ($M = Li, Na, K, \text{ or } Ba$) catalysts for asymmetric Michael reactions. Mechanistic studies on $AlLi-(R)$ -BINOL complex revealed that this complex also acts as a multifunctional heterobimetallic asymmetric catalyst. Furthermore, three-component coupling reactions were achieved by trapping the aluminum enolate intermediate with an aldehyde for the first time. These results will provide a very useful method for the catalytic asymmetric synthesis of complex molecules.

Experimental Procedure

I: To a solution of LiAlH_4 (94.9 mg, 2.5 mmol) in THF (5.0 mL) was added a solution of (*R*)-BINOL (20 mL, 5.0 mmol, 0.25 M in THF) at 0°C. After 30 min stirring at 0°C, this solution of **I** (0.1 M in THF) was directly used for Michael reactions.

7: To a stirred solution of the catalyst **I** (0.05 mmol) in THF (0.5 mL) was successively added cyclohexenone (**6**) (48 μL , 0.5 mmol) and dibenzyl malonate (**4**) (125 μL , 0.5 mmol) at 0°C. After 72 h stirring at room temperature, the reaction mixture was treated with 1 N HCl (2.0 mL) and extracted with EtOAc (3 \times 10 mL). The combined organic extracts were washed with brine, dried (Na_2SO_4), and concentrated to give an oily residue. Purification by flash chromatography (SiO_2 , 25% acetone/hexane) gave the Michael adduct **7** (172.0 mg, 88% yield) in 99% ee. The spectral and analytical data of **7** were in agreement with those of Ref. [2].

13: To a solution of cyclopentenone (**1**) (84 μL , 1.0 mmol), methyl diethylmalonate (**2**) (138 μL , 1.0 mmol) and 3-phenylpropanal (**12**) (158 μL , 1.2 mmol) in THF (1.0 mL) was added a solution of the AlLi -(*R*)-BINOL complex **I** (1.0 mL, 0.1 M in THF, 0.1 mmol) at 0°C. After 36 h stirring at room temperature, the reaction was quenched with 1 N HCl at 0°C and worked up as described for **7**. The oily residue was purified by flash chromatography (SiO_2 , acetone:hexane 1:9) to give **13** (250 mg, 0.64 mmol, 64% yield). The optical purity of **13** was determined by HPLC analysis on a chiral phase (DAICEL CHIRALPAK AD, *i*-PrOH:hexane 1:9). $[\alpha]_D^{25} = +18.3$ ($c = 0.66$ in chloroform) (91% ee); $^1\text{H NMR}$ (270 MHz, CDCl_3): $\delta = 7.40$ –7.13 (m, 5 H), 4.13 (q, $J = 7.0$ Hz, 2 H), 4.11 (q, $J = 7.0$ Hz, 1 H), 4.06 (q, $J = 7.0$ Hz, 1 H), 3.75 (dt, $J = 3.8, 9.6$ Hz, 1 H), 2.97 (dt, $J = 6.3, 8.0$ Hz, 1 H), 2.89–2.76 (m, 1 H), 2.70–2.57 (m, 2 H), 2.35 (dd, $J = 3.8, 6.3$ Hz, 1 H), 2.24 (t, $J = 8.4$ Hz, 2 H), 2.19–2.01 (m, 2 H), 1.87–1.68 (m, 2 H), 1.40 (s, 3 H), 1.22 (t, $J = 7.0$ Hz, 3 H), 1.19 (t, $J = 7.0$ Hz, 3 H); $^{13}\text{C NMR}$ (CDCl_3): $\delta = 218.9, 171.8, 171.6, 141.8, 128.4, 125.9, 72.0, 61.6, 56.9, 55.9, 43.3, 38.5, 36.1, 32.5, 22.9, 18.6, 14.0, 13.9$; IR (neat): $\tilde{\nu} = 3518, 3085, 1728, 1253$ cm^{-1} ; MS: m/z (%) 391 [M^+], 373 [$M^+ - \text{H}_2\text{O}$], 175 (100); anal. calcd for $\text{C}_{22}\text{H}_{30}\text{O}_6$: C, 67.67; H, 7.74; found: C, 67.41; H, 7.70.

16: To a solution of cyclopentenone (**1**) (84 μL , 1.0 mmol), methyl diethylmalonate (**2**) (138 μL , 1.0 mmol) and benzaldehyde (**14**) (122 μL , 1.2 mmol) in THF (1.0 mL) was added a solution of the AlLi -(*R*)-BINOL complex (1.0 mL, 0.1 M in THF, 0.1 mmol) at 0°C. After 72 h stirring at room temperature, the reaction was quenched with 1 N HCl at 0°C and worked up as described for **7**. The crude product was purified by flash chromatography (SiO_2 , acetone:hexane 1:9) to give the diastereomixture of **15** (297 mg, 82%). To a solution of **15** (100 mg, 0.277 mmol) was added molecular sieves (4Å, 200 mg) and pyridinium chlorochromate (PCC, 65 mg, 0.3 mmol) at 0°C. After 1 h of stirring at room temperature, the reaction mixture was diluted with Et_2O (20 mL). The insoluble materials were filtered off through a Celite pad. The filtrate was evaporated to dryness in vacuo. The residue was purified by silica gel chromatography (SiO_2 , acetone:hexane 1:9) to give **16** (100 mg, 100% yield). The optical purity of **16** was determined by HPLC analysis on a chiral phase (DAICEL CHIRALCEL OJ, *i*-PrOH:hexane 1:4). $[\alpha]_D^{25} = -22.0$ ($c = 1.02$ in chloroform) (89% ee); $^1\text{H NMR}$ (270 MHz, CDCl_3): $\delta = 8.00$ –8.05 (m, 2 H), 7.62–7.55 (m, 1 H), 7.52–7.45 (m, 2 H), 4.77 (d, $J = 8.6$ Hz, 1 H), 4.13 (q, $J = 7.3$ Hz, 2 H), 4.03 (dq, $J = 7.3, 10.9$ Hz, 1 H), 3.87 (dq, $J = 7.3, 10.9$ Hz, 1 H), 3.68–3.55 (m, 1 H), 2.48–2.26 (m, 3 H), 2.03–1.03 (m, 1 H), 1.49 (s, 3 H), 1.21 (t, $J = 7.3$ Hz, 3 H), 1.04 (t, $J = 7.3$ Hz, 3 H); $^{13}\text{C NMR}$ (CDCl_3): $\delta = 211.0, 195.1, 171.1, 136.5, 133.3, 129.4, 128.5, 61.4, 59.4, 55.6, 44.4, 38.9, 22.7, 19.7, 13.9, 13.6$; IR (neat): $\tilde{\nu} = 3463, 1728, 1677, 1260$ cm^{-1} ; MS: m/z (%) 361 [$M^+ + 1$], 360 [M^+], 105 (100); anal. calcd for $\text{C}_{20}\text{H}_{24}\text{O}_6$: C, 66.65; H, 6.71; found: C, 66.38; H, 6.76.

Received: August 18, 1995 [Z 83241E]
German version: *Angew. Chem.* 1996, 108, 103–105

Keywords: aldol reactions · aluminum compounds · cascade reactions · catalysis · Michael additions

- [1] a) H. Sasai, T. Suzuki, S. Arai, T. Arai, M. Shibasaki, *J. Am. Chem. Soc.* 1992, 114, 4418; b) H. Sasai, T. Suzuki, N. Itoh, M. Shibasaki, *Tetrahedron Lett.* 1993, 34, 851; c) H. Sasai, N. Itoh, T. Suzuki, M. Shibasaki, *ibid.* 1993, 34, 855; d) H. Sasai, T. Suzuki, N. Itoh, S. Arai, M. Shibasaki, *ibid.* 1993, 34, 2657; e) H. Sasai, W.-S. Kim, T. Suzuki, M. Shibasaki, M. Mitsuda, J. Hasegawa, T. Ohashi, *ibid.* 1994, 35, 6123; f) H. Sasai, T. Suzuki, N. Itoh, K. Tanaka, T. Date, K. Okamura, M. Shibasaki, *J. Am. Chem. Soc.* 1993, 115, 10372; g) H. Sasai, T. Tokunaga, S. Watanabe, T. Suzuki, N. Itoh, M. Shibasaki, *J. Org. Chem.* 1995, 60, 7388; h) H. Sasai, S. Arai, Y. Tahara, M. Shibasaki, *ibid.* 1995, 60, 6656.
- [2] H. Sasai, T. Arai, Y. Satow, K. N. Houk, M. Shibasaki, *J. Am. Chem. Soc.* 1995, 117, 6194, and references therein.
- [3] For the previously known asymmetric Al -BINOL reagents and the related compounds, see: a) R. Noyori, *Tetrahedron* 1994, 50, 4259, and references therein; b) K. Maruoka, H. Yamamoto in *Catalytic Asymmetric Synthesis* (Ed.: I. Ojima), VCH, New York, 1993, p. 413, and references therein.
- [4] Calculated on the assumption of quantitative formation of the catalyst.
- [5] For information on the absolute configurations and enantiomeric excesses of all Michael adducts, see Refs. [2, 6].

- [6] H. Sasai, T. Arai, M. Shibasaki, *J. Am. Chem. Soc.* 1994, 116, 1571, and references therein.
- [7] Prepared from diisobutylaluminum hydride, two equivalents of (*R*)-BINOL, and one equivalent of $\text{NaO}t\text{Bu}$, $\text{KN}(\text{TMS})_2$, and/or $\text{Ba}(\text{O}t\text{Bu})_2$.
- [8] The crystalline **I** (10 mol%) efficiently promoted the Michael reaction of **6** with **4** to give **7** (81% ee, 75% yield).
- [9] The crystalline **II** (10 mol%) efficiently promoted Michael reaction of **6** with **4** to give **7** (89% ee, 67% yield).
- [10] Crystal structure data for **II** ($\text{C}_{40}\text{H}_{24}\text{O}_4\text{AlLi}/\text{C}_6\text{H}_8\text{O}/(\text{C}_4\text{H}_8\text{O})_3$): $T = 223$ K; space group $P1$, $a = 9.510(1)$, $b = 11.063(1)$, $c = 13.052(3)$ Å; $\alpha = 96.42(1)$, $\beta = 95.71(1)$, $\gamma = 64.97(1)^\circ$, $Z = 1$. The structure was solved by direct methods and refined to $R(F) = 0.0520$, $R_w(F) = 0.0490$; further details of the crystal structure investigation are available on request from the Director of the Cambridge Crystallographic Data Centre, 12 Union Road, GB-Cambridge CB2 1EZ (UK) on quoting the full journal citation.
- [11] C. Lambert, P. von R. Schleyer, *Angew. Chem.* 1994, 106, 1187; *Angew. Chem. Int. Ed. Engl.* 1994, 33, 1129.
- [12] The chemical shift of the α -proton of **I** in the ^1H and ^{13}C NMR spectra does not change on mixing **I** and the AlLi -BINOL complex.
- [13] a) D. Canet, J. J. Delpuech, M. R. Khaddar, P. R. Rubini, *J. Mag. Res.* 1973, 9, 329; b) J. J. Delpuech, M. R. Khaddar, A. A. Peguy, P. R. Rubini, *J. Chem. Soc. Chem. Commun.* 1974, 154; c) *J. Am. Chem. Soc.* 1975, 97, 3373.
- [14] The mixture *rac*-**3** does not react with **12** in the presence of 10 mol% of the AlLi -(*R*)-BINOL complex.

Enzymatic Synthesis of Peptides and *Ras* Lipopeptides Employing Choline Ester as a Solubilizing, Protecting, and Activating Group**

Michael Schelhaas, Simone Glomsda, Marion Hänslér, Hans-Dieter Jakubke,* and Herbert Waldmann*

For the synthesis of peptides and sensitive complex peptide conjugates, enzymatic transformations offer advantageous alternatives to classical chemical methods.^{11, 21} Thus, for instance, biocatalyzed deblocking reactions¹¹ have proven their efficiency in the construction of glyco-¹³ and lipopeptides,¹⁴ and protease-mediated peptide synthesis in ice¹⁵ can be used to couple amino acids and peptides without racemization.

However, enzymatic reactions often are hampered or even prevented by the low solubility of the protected peptides in aqueous reaction media, resulting in limited accessibility of the substrates to the biocatalysts. This difficulty can be overcome by employing solubilizing blocking groups.¹⁶ Their removal, however, often requires reaction conditions that cannot be tolerated in the construction of sensitive peptide conjugates. We now report that choline esters are advantageous solubility-enhancing protecting groups for the synthesis of peptides and lipopeptides under extremely mild conditions.

Peptide choline esters (peptide-Cho) **1** (Scheme 1) can be constructed by treating the corresponding 2-bromoethyl esters with trimethylamine.¹⁷ Alternatively, they are accessible by condensation of amino acid choline ester hydrobromides with C-terminally deblocked peptides (vide infra; Scheme 2).

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**] This work was supported by the Deutsche Forschungsgemeinschaft, the Bundesministerium für Forschung und Technologie, and the Fonds der Chemischen Industrie. We are grateful to Degussa AG for the donation of chemicals. We thank Mrs. R. Schauf and Mrs. H. Späth for assistance with the experiments.

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